

## DNAfectin™ Plus Transfection Reagent

Cat. No.	Description	Quantity
G2500	DNAfectin™ Plus	1.0 ml



### Storage Conditions

- Store at 4°C. Do not freeze



### Required Materials

- DNA (0.2-16 µg)
- Serum-free, antibiotic-free medium
- Microcentrifuge tubes



### Timing

Preparation: 10 minutes  
 Incubation: 20 minutes  
 Total Incubation: 12-16 hours



### Description

abm's DNAfectin™ Plus is a nanoparticle-based, nonliposomal formulation that enables the efficient transfection of plasmid DNA and short oligonucleotides into a broad range of cells with minimal cytotoxicity. This simple protocol does not require the removal of serum or culture medium, resulting in less variability and low risk of contamination. DNAfectin™ Plus has been shown to transfect a wide variety of primary, adherent and suspension cell lines with high efficiency.



### Transfection Optimization

To achieve the maximum transfection efficiency and low cytotoxicity, optimize the transfection conditions by varying cell density along with DNA and DNAfectin™ Plus concentrations. Optimal results have been observed when cells are 80-90% confluent and DNA(µg): DNAfectin™ Plus (µl) ratios are 1:1 to 1:5.

**Table 1: Reagent Quantities for Different Culture Vessels**


Culture Vessel	Volume of plating medium per well	DNA(µg)	DNAfectin™ Plus (µl)	Transfection medium volume
24-well	500µl	0.2-0.4µg	0.6-1.2µl	50µl
12-well	1ml	0.5-0.8µg	1.5-2.5µl	100µl
6-well	2ml	1.0-2.0µg	3-6µl	200µl
35mm	2ml	1.0-2.0µg	3-6µl	200µl
60mm	5ml	3.0-6.0µg	10-20µl	300µl
10cm	10ml	8.0-16.0µg	25-50µl	500µl

## Transfection Protocol

Use the following conditions as guidelines to transfect mammalian cells in a 6-well or 35mm dish format. For other culture vessels, please refer to Table 1.

1. **Plating Cells:** 18 to 24 hours prior to transfection, seed the cells at a density such that they are in optimal culture conditions. Incubate the cells at 37°C in a CO<sub>2</sub> incubator until the cells are 70% to 90% confluent at the time of transfection.
2. For each transfection sample, prepare the DNAfectin™ Plus-DNA complexes as follows:
  - a) Add 2.0 µg of DNA into 200 µl of serum-free, antibiotic-free medium.
  - b) Warm the DNAfectin™ Plus to room temperature and vortex gently before use.
  - c) Add 6.0 µl of the DNAfectin™ Plus into the DNA solution from step a). Pipette up and down gently several times to mix the solution completely.
  - d) Incubate for 20 minutes at room temperature to form the DNAfectin™ Plus-DNA complexes. Complexes are stable at room temperature for 3-5 hours.
3. Transfer the DNAfectin™ Plus-DNA solution to the cultured cells drop-by-drop to different areas of the culture dish. Gently rock the culture vessel back-and-forth and side-to-side to evenly distribute the complexes.
4. Incubate for 12-16 hours. It is not necessary to change the culture medium after transfection with DNAfectin™ Plus, however, culture medium may be changed between 6-24 hours after transfection for sensitive cell lines.
5. Monitor transfection efficiency 24-72 hours post-transfection using relevant assays.

**Figure 1: Performance Data - DNAfectin™ Plus Transfection Reagent**



### Notices and Disclaimers

abm products are intended for laboratory research purposes only, unless noted otherwise. They are not intended for use in humans

