

Product Stability Report

Blastaq™ 2X qPCR MasterMix

Cat. No. G891, G892



Applied Biological Materials Inc.
3671 Viking Way, Unit 1
Richmond, BC, Canada
V6V 2J5
www.abmgood.com

Blastaq™ 2X qPCR MasterMix Product Stability Report

Product Stability Report

Purpose

To test the stability of **Blastaq™ 2X qPCR MasterMix** under several conditions

Method

Blastaq™ 2X qPCR MasterMix was subjected to the following conditions:

- Freeze Thaw (20 cycles)
- 4°C (1 Day, 3 Days, 7 Days)
- 20°C (1 Day, 3 Days, 7 Days)
- 55°C (1 Day, 3 Days, 7 Days)

The qPCR reactions were assembled according to G891 Product Datasheet. Negative control reactions (NTCs) received an equal volume of nuclease-free water instead of the template. The freeze-thaw results were compared to an identical tester that did not undergo any treatment. The results from the “Elevated Temperature” Test were compared against **Blastaq™ 2X qPCR MasterMix** stored under standard conditions (-20 °C). Each dilution was performed in duplicate. After data collection, the melting curve peaks were recorded and analyzed, and any changes in Ct values among the samples were examined.

Result

Blastaq™ 2X qPCR MasterMix samples with comparable qPCR amplification activity to the control are noted as a “Pass” and any condition that is showing less activity is noted with “Decrease in activity”.

| Sample Name | 20 X | 1 Day | 3 Days | 7 Days |
|-------------|------|-------|--------|-------------|
| Freeze-thaw | Pass | -- | -- | -- |
| 4 °C | -- | Pass* | Pass* | Pass* |
| 20 °C | -- | Pass* | Pass* | Pass* |
| 55 °C | -- | Pass* | Pass* | No activity |

* Decrease in fluorescence may be observed

Conclusion

Blastaq™ 2X qPCR MasterMix remains stable and retains full functionality after 7 days at both 4°C and 20°C, as well as after 20 cycles of freeze-thaw. However, repeated freeze-thaw cycles and storage at elevated temperatures are not recommended, as they may result in a reduced fluorescence signal. Furthermore, incubation at 55°C should not exceed 3 days, as this will result in complete loss of activity.

The product has been tested and shown to be in compliance with all specifications.