



Safe-Red™ Nucleic Acid Stain

Cat. No. G108-R

Store at 18-25°C.

Product Description

Safe-Red™ is a new and safe class of gel loading dye for the visualization of double-stranded DNA, single-stranded DNA, and RNA in agarose gels. The dye was developed to replace Ethidium Bromide (EtBr, a potent mutagen), commonly used in gel electrophoresis for visualization of nucleic acids in agarose gels. Safe-Red™ is non-carcinogenic by the Ames-test. The results are negative in both the mouse marrow chromophilous erythrocyte micronucleus and mouse spermary spermatocyte chromosomal aberration tests. Safe-Red™ loading dye does not require the addition of any stain to the molten agarose, simplifying workflow processes. The product is provided in a ready to use 6X format and is added directly to samples prior to loading onto the gel. Safe-Red™ loading dye eliminates contamination risk of glassware or agarose gel apparatus as associated with EtBr.

Cat. No.	Product	Concentration	Format	Quantity
G108-R	Safe-Red™	6X	Loading Dye	1.0 ml

Protocol

1. Prepare agarose gel solution, mix gently to avoid air bubbles.
2. Let the solution cool to 60-70°C and cast the gel.
3. Mix samples and DNA marker with Safe-Red™ loading dye at a 1:5 (dye : sample) dilution ratio, and load onto gel.
4. Following electrophoresis, view the results under UV or blue LED light.