


# Product Testing Report

---



## RNaseOFF: Product Efficiency


### Purpose

To compare the performance of RNaseOFF Ribonuclease Inhibitor to two other leading brand competitor products.

### Method

To compare the performance, competitor products (Competitor I and Competitor P) were purchased and quantified using SDS-PAGE. Next, all three proteins were normalized to 0.15 mg/ml and a premix of substrate (HEK293T RNA) (1 µg) and ribonuclease inhibitor (1.5 µg) was assembled and aliquoted into PCR tubes, followed by the addition of RNase A (AG Scientific) at increasing concentrations (0-10 ng). Positive (10 ng RNase A added) and negative controls (no RNase A added) demonstrated how the substrate would look if no ribonuclease inhibitor was added. The reactions were incubated (15 min, 37 °C) and subjected to agarose gel electrophoresis (10 µl load).

**Figure 1 – RNaseOFF offers higher substrate protection compared to two leading competitor products when 10 ng of RNase A is added.** Intact RNA should exhibit two clear bands, the 18S and 28S ribosomal RNA band, on the agarose gel. Degraded RNA appears as low molecular weight smear.



# Product Testing Report

---

## Result

Equal performance (intact HEK293T RNA) was observed between all three products in the presence of 0-5 ng RNase A. However, RNaseOFF Ribonuclease Inhibitor was superior when 10 ng of RNase A was added: highlighted lanes demonstrate significant substrate digestion in case of the two competitor products but almost complete integrity in case of RNaseOFF.

## Conclusion

RNaseOFF Ribonuclease Inhibitor exhibits better activity compared to Competitor I and Competitor P at an identical concentration.

## RNaseOFF: Product Stability

### Purpose


To assess the stability of RNaseOFF Ribonuclease Inhibitor at room temperature (RT) 20 °C.

### Method

RNaseOFF Ribonuclease Inhibitor (1.5 µg) was incubated at 20 °C and 4 °C (to mimic transportation conditions) for 2, 4, and 7 days. Collected testers were frozen, subsequently thawed all at once, and added to the substrate (HEK293T RNA) (1 µg). The premixes were subjected to incubation with RNase A (AG Scientific, 5 ng) and agarose gel electrophoresis (10 µl load).

Figure 2 – RNaseOFF retains its stability at 4 degrees and 20 degrees (RT) for up to 7 days. It successfully inhibits the RNase A in all samples tested

T (°C)	-20°C	-20°C	RT	4°C	RT	4°C	RT	4°C	-20°C	-20°C
Day	-	-	2	2	4	4	7	7	-	-
RNase A (ng)	0 ng	5 ng	5 ng	5 ng	5 ng	5 ng	5 ng	5 ng	0 ng	0 ng



abm RNaseOFF

No abm  
RNaseOFF

## Result

No RNA degradation was observed if the RNaseOFF remained at 20 °C or 4 °C for up to 7 days. The product maintained its stability and protected the substrate in the presence of 5 ng of RNase A.

**Note that RNase A is a highly stable enzyme that remains functional when heated to 100 °C for 15 minutes** (Narahashi, T., & Conn, P. M. Ion channels of excitable cells. (2013). Academic Press; RNase A: Frequently Asked Questions. (2019, May). AG Scientific).

Therefore, it had been actively inhibited by the RNase OFF Ribonuclease Inhibitor for the whole duration of the experiment.

## Conclusion

While we recommend storing RNase OFF Ribonuclease Inhibitor at -20 °C, the product may be used for experiments that require prolonged substrate incubations at 2-6 °C.