



## BlasTaq™ 2X qPCR MasterMix

### Cat. No. G891, G892

Store at -20°C.

### Product Description

**BlasTaq™ 2X qPCR MasterMix** provides a **dye-based** real-time qPCR analysis of DNA samples. This ready-to-use qPCR MasterMix contains **abm's** strategically-engineered, next generation Taq Polymerase, BlasTaq™, leading to **rapid extension rates and robust performance**. With specialized reaction conditions, this polymerase provides increased processivity, yields, and sensitivity, while shortening reaction times by up to 70%, compared to wild-type Taq DNA polymerase. BlasTaq™ has 5'-3' polymerase and 5'-3' exonuclease activities, lacks 3'-5' exonuclease activity, and produces 3'-dA-tailed amplicons. qPCR products made with BlasTaq™ can be used with TA cloning vectors. MasterMix contains dye comparable to SYBR Green™ and EvaGreen™. ROX Reference Dye is provided separate from the MasterMix, making this kit universally compatible with most qPCR instruments.

Cat. No.	Product Component	Quantity	Part No.
G891	BlasTaq™ 2X qPCR MasterMix	500 rxn (4 x 1.25 ml)	G891-1
	ROX Reference Dye	50 µl	P102
G892	BlasTaq™ 2X qPCR MasterMix	2,500 rxn (25 ml)	G892-1
	ROX Reference Dye	240 µl	P103

### Protocol

MasterMix contains dye comparable to SYBR Green™ and EvaGreen™. ROX Reference Dye is provided separate from the MasterMix, making this kit universally compatible with most qPCR instruments.

See **Rox Machine Compatibility** on our product page under the Documents tab on our website.

The recommended amount of ROX Reference Dye to be added into the MasterMix may vary depending on the qPCR machine type:

- No ROX equipment: Not needed.
- Low ROX equipment: 1 µl/1.25 ml or 22.5 µl/25 ml MasterMix.
- High ROX equipment: 11.5 µl/1.25 ml or 225 µl/25 ml MasterMix

1. Mix individual components before use and assemble reaction on ice.

Product Component	Volume
BlasTaq™ 2X qPCR MM <sup>1</sup>	10 µl
Forward Primer (10 µM)	0.5 µl
Reverse Primer (10 µM)	0.5 µl
Template DNA	Variable (100 ng genomic DNA)
Nuclease-free H <sub>2</sub> O	Up to 20 µl

<sup>1</sup> The reaction buffer contains 1.5 mM Mg<sup>2+</sup>

2. Gently mix the reaction components, and briefly centrifuge. Use thermocycling conditions below.

Step	Temperature	Duration		Cycle(s)
		Standard	Fast	
Enzyme Activation	95°C	3 min	3 min	1
Denaturation	95°C	15 sec	1 sec	40
Annealing/Extension	60°C	1 min	10 sec	
Melting Curve	Refer to specific guidelines for instrument used			

### General Notes

- Specialized buffer for higher yields, sensitivity, and specificity compared to wild-type Taq polymerase.
- Ideally start the qPCR as soon as the reaction mixture is prepared. If not possible, keep the reaction mixture on ice until starting the qPCR.
- Use the standard thermocycling condition with miRNA cDNA templates or any other appropriate applications.