



## ViralEntry™ Transduction Enhancer

### Cat. No. G515

Store at 4°C.

### Product Description

As an industry leader in recombinant viral vectors, abm has 15+ years of experience in how to enhance virus vector transduction efficiency. Over the years, we have developed multiple formulations of viral enhancing agents useful for this purpose. A remarkable breakthrough by our top scientists in 2020 is the development of the most efficient viral enhancing agent, **ViralEntry™, that is not only effective in enhancing lentiviral vector transduction efficiency significantly, but also increases transduction efficiency in diverse cell types** including primary T lymphocytes by over 10X. This is superior to any similar product on the market.

With this incredible reagent, **you will be able to transduce any cell types you have.** We guarantee that this is the most effective viral transduction reagent in the industry and will gladly refund you if you find better transduction enhancers than our ViralEntry™ formulation. Use it and you will enjoy it

Product Component	Quantity
ViralEntry™	1.0 ml

### Protocol

#### A) For adherent and/or suspension cells (direct culture):

Note: Protocol below is for a 6-well plate.

- Day 1.** Plate your cells of interest in a 6-well plate 24 hours before infection with a density of  $2 \times 10^5$  cells/well.
- Day 2.** Add ViralEntry™ at 1:100 to culture medium, and mix gently. Infect each well with recombinant virus at optimal MOI. Incubate cells at appropriate temperature (°C) and CO<sub>2</sub> (%) conditions.  
$$\text{MOI} = \frac{\text{Infection Units of virus used}}{\text{number of cells}}$$
- Day 4.** Change growth medium. If the recombinant virus contains a drug resistance gene for stable cell clone generation, subculture cells to 10 cm dish for drug selection. If the recombinant virus contains a fluorescent tag, you can evaluate transduction efficiency by checking signals under the fluorescence microscope.

#### B) For suspension cells (spinoculation):

Note: Protocol below is for a 24-well plate.

##### 1. Day 1.

- Dilute ViralEntry™ 1:100 with appropriate growth medium. Mix gently.
- Add 400 µl of 1:100 enhancer mix to one well of a 24 well plate, and place in CO<sub>2</sub> incubator.
- Add  $1 \times 10^5$  viable cells in a sterile 1.5 ml Eppendorf tube, spin down at 300 xg for 5 minutes, aspirate supernatant, and resuspend pellet with 100 µl of the 1:100 enhancer mix.
- Add recombinant virus at optimal MOI, and mix by gentle tapping.
- Keep at room temperature for 10 minutes while gently tapping the tube every 2-3 minutes.
- Spin down at 1700 xg for 20 minutes.
- Gently resuspend the cell pellet and transfer suspension to the 24-well plate containing 400 µl diluted ViralEntry™ (step b). Incubate the cells at appropriate temperature (°C) and CO<sub>2</sub> (%) conditions.

- Day 3.** Change growth medium. If the recombinant virus contains a drug resistance gene, begin drug selection by replacing media with drug containing media every 3-4 days until resistant cells can be identified. If the recombinant virus contains a fluorescent tag, you can evaluate transduction efficiency by checking signals under the fluorescence microscope.

### General Notes

- Spinoculation protocol is recommended for T cells, and other difficult to transduce suspension cells.