

High Speed Micro-centrifuge with 24 x 1.5 ml Rotor

Cat. No. Q5140

Buy now for only **\$1500** USD

Powerful high speed micro-centrifuge at an incredible price

- Comes with a fixed angle rotor that can hold 24 x 1.5 ml tubes
- Sturdy, stainless steel chamber provides reliable performance
- Compact and portable design saves you valuable lab space

Rapid acceleration with exceptional control

- Reach up to 20,000 x g in just 16 seconds
- Speed range of 100-15,000 rpm, adjustable in increments of 100 rpm

Simplify routine operations with easy-to-use controls

- Turn a knob to quickly change parameters
- Push a button to easily convert between RPM and RCF



Goes up to 15,000 rpm
or 20,375 x g



Holds 24 x 1.5 ml
tubes



1 Year Manufacturer's
Warranty

SPECIFICATIONS

Max Capacity..... 24 x 1.5 ml tubes
 Max Speed 15,000 rpm
 Max RCF..... 20,375 x g
 Max Accel/Decel Time 16 s / 18 s
 Time Range 1-99 minutes
 Noise Level..... ≤60 dB(A)

Power Supply 120 V
 Dimensions (D x W x H) 280 x 360 x 250 mm
 Total Weight..... 23 kg
 Auto Lid Release Yes, with manual release cord
 Safety Alarms Automatic imbalance detection
 Certifications..... ISO 13485:2003, CE, REACH, RoHS