



Product Datasheet

Cat. No.
TM040

Name
Keratinocyte Medium Kit

Description

abm's Keratinocyte Medium Kit is a low serum medium developed for the cultivation of keratinocytes in culture. The media contains proprietary formulation of growth factors and supplements to support optimal keratinocyte growth. This media supports robust growth of keratinocytes while maintaining characteristic morphology and markers of keratinocytes in culture.

Pack Size

The kit has following components:

- Basal Media (500 ml)
 - Keratinocyte Media Supplements
-

Product use

To make complete medium, add the entire content of the Keratinocyte Media Kit supplements, 2% FBS and appropriate antibiotics into the Basal Media. Each lot has been fully tested and optimized for the ability to support the growth of the intended primary cells, in addition to being QC tested for sterility, pH, osmolality and endotoxin levels.

Storage

The shelf life of Basal Media is stored at 4-8°C while the Keratinocyte Media Kit supplement is stored at -20°C. The media is good to use until listed expiry date on bottle. Media expires 30 days after supplementation.

Notices and Disclaimers

1. **abm** products are intended for laboratory research purposes only, unless noted otherwise. They are not intended for use in humans.
 2. **abm** makes no warranties or representations as to the accuracy of the information on this site. Citations from literature and provided for informational purposes only. **abm** does not warrant that such information has been shown to be accurate.
-

Caution: This product is for research use only and is not intended for therapeutic or diagnostic applications.
Please contact a technical service representative for more information (1-866-757-2414).

1 - 3 6 7 1 Viking Way
Richmond BC, Canada V6V 2J5
Tel : 6 0 4 - 2 4 7 - 2 4 1 6
Fax : 6 0 4 - 2 4 7 - 2 4 1 4
www.abmGood.com