



Applied Biological Materials Inc

Tel: 1-866-757-2414
Email: info@abmGood.com
Website: www.abmGood.com

Sodium Pyruvate Solution

Store at 2-8°C

Cat. No.	Description	Quantity
TM057	Sodium Pyruvate Solution	100 ml

Product Description

Sodium pyruvate is used in cell culture as a source of carbon. It is also involved in amino acid metabolism and initiates the Krebs cycle. Since cells can produce sodium pyruvate as an intermediate product during glycolysis, it is a required supplement to all cell culture.

Usage

At concentration of 100 mM, **abm's** Sodium Pyruvate Solution is ready to be used for your cell culture needs. Usually used at final concentration of either 1 mM or 2 mM in medium.

Storage Conditions

Store at 2-8°C.

Shelf Life

18 months from the date of shipping, if stored and handled correctly. Do not freeze.

For laboratory research only. Not for clinical applications.
For technical questions, phone the ABM helpline at 1-866-757-2414

or visit our website at www.abmGood.com